

**DESIGNED
TO LAST**



mini 

The Compact Automated Immunoassay Analyzer

- Easy to use
- Cost Effective
- Reliable

With over 40 years of experience in the "in vitro" diagnostics market place, bioMérieux also sets the standard of excellence in fully automated immunoanalysis. With over 18,000 VIDAS® and miniVIDAS installed worldwide, these systems have become a reference in the immunoassay field.

mini **VIDAS**
DESIGNED TO LAST



Tailor-made to fit your lab

Rapid turnaround times and independent testing design make miniVIDAS the ideal instrument for both routine batches and stat tests.

The two 6-test sections enable you to process different parameters simultaneously with a throughput rate of 80 samples per day.

The all inclusive concept that ensures effective cost control

All required consumables are provided in the reagent kit: calibration, control, substrate, diluents and master curve.

The reagent strip contains all the ready-to-use reagents and each strip is bar-coded for automated identification of the test.

Designed to reach high quality results

The ready-to-use disposable VIDAS test is completely self-contained. There are no sampling needles, syringes or tubing and no contact between reagent and instrument.

This unique design virtually ensures the robustness of the system, minimum maintenance and prevents cross-contamination in particular with analyses which require very low limits of detection.

Dimensions:

Height 17.7"
Width 22.6"
Depth 22"
Weight 82 lbs.

bioMérieux, Inc.
100 Rodolphe Street
Durham, NC 27712
U.S.A.
Tel: (800) 682 2666
Fax: (800) 968 9494

bioMérieux Canada, Inc.
7815 Henri-Bourassa Boulevard West
St. Laurent (Québec)
Canada H4S 1P7
Tel: (1) 514 336 7321
Fax: (1) 514 807 0015

www.biomerieux-usa.com


BIOMÉRIEUX