

Performa®

*Merit's Diagnostic Cardiology Catheters have
an Ideal Combination of Torque, Push & Flow*



It All Adds Up

MERIT MEDICAL DIAGNOSTIC CARDIOLOGY CATHETERS

Introducing Merit Medical's line of Diagnostic Cardiology catheters. Designed with all the torque you need, improved shaft strength for better pushability and large inner lumens for increased flow rates. Now you can reach your goals in today's cardiology environment with exceptional products from Merit Medical.

Performa®

NYLON PEBAX MATERIAL

- Kink Resistant
- Optimal surface smoothness
- Thermal stability at body temperature
- Thin wall design for improved flow rate

RADIOPACITY

- Good shaft visualization under fluoroscopy
- Tip opacity ensures accurate placement and visualization

WINGED POLYCARBONATE HUB

- Enhanced handling and control

WIRE BRAID DESIGN

- Greater torque for maximum control
- Improved longitudinal stability
- High degree of maneuverability
- Provides superior pushability

TIP DESIGNS

- Excellent tip memory
- Reforms accurately time after time

BUMPER TIP

- Added softness to the 4F Performa catheter line for patient safety

FLAT WIRE BRAID

- Increased Torque
- Greater Stability
- Improved Kink Resistance



Proven Performer | Great Radiopacity | Good Torque

4F PERFORMA CATHETER LINE

Accepts .038" (0.97mm) guide wire for physician convenience

Bumper tip softness added for safe placement of tip into ostium

Large inner diameter .042" (1.07mm) provides excellent flow rates

INNER-LUMEN SIZES

- 4F .042" (1.07mm)
- 5F .046" (1.17mm)
- 5F .052" (1.32mm) high flow
- 6F .054" (1.37mm)
- 6F .059" (1.49mm) high flow
- Yields higher flow rates at lower pressures
- Superior guide wire trackability

STRAIN RELIEF

- Increased strength at catheter/hub interface
- Color coded for easy identification

PRESSURE

- Rated at 1200psi (81.6 bar)

PERFORMA CARDIOLOGY CATHETERS

catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--


Judkins Left

7700-10	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL3.5	bumper tip
7700-20	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL4.0	bumper tip
7700-30	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL5.0	bumper tip
7700-40	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL6.0	bumper tip
7521-13	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL3.5	
7521-21	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL4.0	
7521-F1	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL4.5	
7521-33	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL5.0	
7521-43	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL6.0	
7521-53	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL3.5	
7521-63	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL4.0	
7521-73	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL5.0	
7521-83	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL6.0	
7601-13	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL3.5-High Flow	non-braided
7601-21	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL4.0-High Flow	non-braided
7601-33	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL5.0-High Flow	non-braided
7601-43	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL6.0-High Flow	non-braided
7501-13	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL3.5	
7501-21	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL4.0	
7501-H1	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL4.5	
7501-33	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL5.0	
7501-43	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL6.0	
7501-53	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL3.5	
7501-63	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL4.0	
7501-73	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL5.0	
7501-83	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL6.0	
7574-13	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL3.5-High Flow	non-braided
7574-21	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL4.0-High Flow	non-braided
7574-33	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL5.0-High Flow	non-braided
7574-43	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL6.0-High Flow	non-braided



Judkins Left Short-tip


7501-A3	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL3.5 ST	
7501-B1	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL4.0 ST	
7501-C3	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL5.0 ST	
7501-D3	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL6.0 ST	



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--


Judkins Right

7701-10	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR3.5	bumper tip
7701-20	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR4.0	bumper tip
7701-30	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR5.0	bumper tip
7701-40	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR6.0	bumper tip
7523-13	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR3.5	
7523-21	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR4.0	
7523-33	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR5.0	
7523-43	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR6.0	
7523-53	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR3.5	
7523-63	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR4.0	
7523-73	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR5.0	
7523-83	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR6.0	
7503-13	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR3.5	
7503-21	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR4.0	
7503-33	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR5.0	
7503-43	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR6.0	
7503-53	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR3.5	
7503-63	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR4.0	
7503-73	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR5.0	
7503-83	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR6.0	




Judkins Right Short-tip

7523-E3	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR3.5 ST	
7523-F1	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR4.0 ST	
7523-G3	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR5.0 ST	
7523-H3	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR6.0 ST	
7503-A3	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR3.5 ST	
7503-B1	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR4.0 ST	
7503-C3	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR5.0 ST	
7503-D3	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR6.0 ST	



Judkins Right Modified


7502-13	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR3.5	MOD
7502-23	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR4.0	MOD
7502-33	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR5.0	MOD
7502-43	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR6.0	MOD



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--


Amplatz Left

7703-10	4	.042"/1.07mm	100cm	5.0cm	.038" (0.97mm)	AL1	bumper tip
7703-20	4	.042"/1.07mm	100cm	2.2cm	.038" (0.97mm)	AL2	bumper tip
7703-30	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	AL3	bumper tip
7703-40	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	AL4	bumper tip
7527-13	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AL1	
7527-23	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AL2	
7527-33	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	AL3	
7527-43	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	AL4	
7507-13	6	.054"/1.37mm	100cm	5.0cm	.038" (0.97mm)	AL1	
7507-23	6	.054"/1.37mm	100cm	5.0cm	.038" (0.97mm)	AL2	
7507-33	6	.054"/1.37mm	100cm	2.4cm	.038" (0.97mm)	AL3	
7507-43	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	AL4	



Amplatz Right

7702-40	4	.042"/1.07mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	bumper tip
7702-10	4	.042"/1.07mm	100cm	3.3cm	.038" (0.97mm)	AR1	bumper tip
7702-20	4	.042"/1.07mm	100cm	5.0cm	.038" (0.97mm)	AR2	bumper tip
7702-30	4	.042"/1.07mm	100cm	5.0cm	.038" (0.97mm)	AR3	bumper tip
7526-11	5	.046"/1.17mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	
7528-13	5	.046"/1.17mm	100cm	3.3cm	.038" (0.97mm)	AR1	
7528-23	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AR2	
7528-33	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AR3	
7506-11	6	.054"/1.37mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	
7508-13	6	.054"/1.37mm	100cm	3.3cm	.038" (0.97mm)	AR1	
7508-23	6	.054"/1.37mm	100cm	3.5cm	.038" (0.97mm)	AR2	
7508-33	6	.054"/1.37mm	100cm	4.0cm	.038" (0.97mm)	AR3	



Left Coronary Bypass

7709-10	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	LCB	bumper tip
7534-13	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	LCB	
7510-13	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	LCB	



Right Coronary Bypass

7709-20	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	RCB	bumper tip
7534-23	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	RCB	
7510-23	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	RCB	



Progressive Right Coronary (NON TORQUE)

7707-10	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	PRC	bumper tip
7600-11	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	PRC	non-braided
7532-13	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	PRC	
7572-11	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	PRC	non-braided
7512-13	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	PRC	
7512-23	6	.054"/1.37mm	125cm	2.5cm	.038" (0.97mm)	PRC	



Internal Mammary

7706-10	4	.042"/1.07mm	100cm	3.3cm	.038" (0.97mm)	IMA	bumper tip
7531-30	5	.046"/1.17mm	65cm	3.3cm	.038" (0.97mm)	IMA	
7531-13	5	.046"/1.17mm	100cm	3.3cm	.038" (0.97mm)	IMA	
7511-13	6	.054"/1.37mm	100cm	3.3cm	.038" (0.97mm)	IMA	



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	# side ports	small pigtail curve	angle length
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--------------	---------------------	--------------

Straight Cardiac Pigtail

7704-70	4	.042"/1.07mm	80 cm	6.5 cm	.038" (0.97mm)	Braided PIG	5		
7704-10	4	.042"/1.07mm	110 cm	6.5 cm	.038" (0.97mm)	Braided PIG	5		
7776-11	5	.046"/1.17mm	90cm	5.0cm	.038" (0.97mm)	Braided PIG	5	yes	
7776-21	5	.046"/1.17mm	110cm	5.0cm	.038" (0.97mm)	Braided PIG	5	yes	
7776-B1	5	.052"/1.32mm	110cm	4.6cm	.038" (0.97mm)	Non-Braided PIG	5	yes	
7602-11	5	.052"/1.32mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided PIG	5		
7529-11	5	.046"/1.17mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	5		
7529-43	5	.046"/1.17mm	125cm	7.5cm	.038" (0.97mm)	Braided PIG	5		
7569-11	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided PIG	5		
7509-11	6	.054"/1.37mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	5		
7777-B1	6	.059"/1.49mm	110cm	4.8cm	.038" (0.97mm)	Non-Braided PIG	5	yes	
7777-21	6	.054"/1.49mm	110cm	4.6cm	.038" (0.97mm)	Braided PIG	5	yes	
7509-43	6	.054"/1.49mm	125cm	7.5cm	.038" (0.97mm)	Braided PIG	5		
7569-43	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided PIG	5		



PIGTAIL

145 Degree Cardiac Pigtail


7704-20	4	.042"/1.07mm	110cm	6.5cm	.038" (0.97mm)	Braided	5		5.5cm
7776-C3	5	.052"/1.32mm	110cm	4.6cm	.038" (0.97mm)	Non-Braided	5	yes	5.5cm
7776-33	5	.046"/1.17mm	110cm	5.0cm	.038" (0.97mm)	Braided	5	yes	5.5cm
7529-53	5	.046"/1.17mm	125cm	7.5cm	.038" (0.97mm)	Braided	5		7.5cm
7602-23	5	.052"/1.32mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided	5		7.5cm
7529-23	5	.046"/1.17mm	110cm	7.5cm	.038" (0.97mm)	Braided	5		7.5cm
7569-23	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided	5		5.5cm
7509-23	6	.054"/1.37mm	110cm	7.5cm	.038" (0.97mm)	Braided	5		5.5cm
7777-C3	6	.059"/1.49mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided	5	yes	5.5cm
7777-33	6	.054"/1.37mm	110cm	4.6cm	.038" (0.97mm)	Braided	5	yes	5.5cm
7509-53	6	.054"/1.37mm	125cm	7.5cm	.038" (0.97mm)	Braided	5		5.5cm
7569-53	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided	5		5.5cm



145° PIGTAIL

155 Degree Cardiac Pigtail

7704-30	4	.042"/1.07mm	110cm	6.5cm	.038" (0.97mm)	Braided	5		5.5cm
7602-33	5	.052"/1.32mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided	5		7.5cm
7529-33	5	.046"/1.17mm	110cm	7.5cm	.038" (0.97mm)	Braided	5		7.5cm
7529-63	5	.046"/1.17mm	125cm	7.5cm	.038" (0.97mm)	Braided	5		7.5cm
7776-D3	5	.052"/1.32mm	110cm	4.6cm	.038" (0.97mm)	Ultra High Flow	5	yes	5.5cm
7776-43	5	.046"/1.17mm	110cm	5.0cm	.038" (0.97mm)	Braided	5	yes	5.5cm
7569-33	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided	5		5.5cm
7509-33	6	.054"/1.37mm	110cm	7.5cm	.038" (0.97mm)	Braided	5		5.5cm
7509-63	6	.054"/1.37mm	125cm	7.5cm	.038" (0.97mm)	Braided	5		5.5cm
7569-63	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided	5		5.5cm
7777-D3	6	.059"/1.49mm	110cm	4.8cm	.038" (0.97mm)	Non-Braided	5	yes	5.5cm
7777-43	6	.054"/1.37mm	110cm	4.6cm	.038" (0.97mm)	Braided	5	yes	5.5cm




155° PIGTAIL

catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	# side ports	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--------------	--

Multipurpose

7708-30	4	.042"/1.07mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2	bumper tip
7708-40	4	.042"/1.07mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2	bumper tip
7708-10	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-	bumper tip
7708-20	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2	bumper tip
7708-50	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-	bumper tip
7708-60	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2	bumper tip
7533-53	5	.046"/1.17mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2	
7533-63	5	.046"/1.17mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2	
7533-13	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-	
7533-23	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2	
7533-33	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-	
7533-43	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2	
7513-A3	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI A1	-	
7513-13	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2	
7513-D3	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI B1	-	
7513-33	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2	
7513-B3	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-	
7513-23	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2	
7513-E3	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-	
7513-43	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2	
7513-C3	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI A1	-	
7513-53	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI A2	2	
7513-F3	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI B1	-	
7513-63	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI B2	2	




Sones

7705-10	4	.042"/1.07mm	80cm	7.3cm	.038" (0.97mm)	SONES	4	bumper tip
7705-20	4	.042"/1.07mm	100cm	7.3cm	.038" (0.97mm)	SONES	4	bumper tip
7530-13	5	.046"/1.17mm	80cm	7.3cm	.038" (0.97mm)	SONES	4	
7530-23	5	.046"/1.17mm	100cm	7.3cm	.038" (0.97mm)	SONES	4	



Ultimate Radial

410038ULT1	4	.042"/1.07mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 1	-	bumper tip
4100381ULT1	4	.042"/1.07mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 1	1	bumper tip
510038ULT1	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 1	-	bumper tip
5100381ULT1	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 1	1	bumper tip
610038ULT1	6	.054"/1.37mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 1	-	bumper tip
6100381ULT1	6	.054"/1.37mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 1	1	bumper tip
410038ULT2	4	.042"/1.07mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 2	-	bumper tip
4100381ULT2	4	.042"/1.07mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 2	1	bumper tip
510038ULT2	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 2	-	bumper tip
5100381ULT2	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 2	1	bumper tip
610038ULT2	6	.054"/1.37mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 2	-	bumper tip
6100381ULT2	6	.054"/1.37mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 2	1	bumper tip
410038ULT3	4	.042"/1.07mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 3	-	bumper tip
4100382ULT3	4	.042"/1.07mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 3	2	bumper tip
510038ULT3	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 3	-	bumper tip
5100382ULT3	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 3	2	bumper tip
610038ULT3	6	.054"/1.37mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 3	-	bumper tip
6100382ULT3	6	.054"/1.37mm	100cm	2.0cm	.038" (0.97mm)	ULTIMATE 3	2	bumper tip



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--

Internal Mammary Left Brachial/Radial

7242-10	4	.042"/1.07mm	100cm	4.2cm	.038" (0.97mm)	IMLBR1	bumper tip
7242-20	4	.042"/1.07mm	100cm	4.2cm	.038" (0.97mm)	IMLBR2	bumper tip

Left Brachial/Radial


7240-10	4	.042"/1.07mm	100cm	3.4cm	.038" (0.97mm)	LBR1	bumper tip
7240-20	4	.042"/1.07mm	100cm	3.4cm	.038" (0.97mm)	LBR2	bumper tip
7240-30	4	.042"/1.07mm	100cm	6.0cm	.038" (0.97mm)	LBR3	bumper tip

Right Brachial/Radial

7241-10	4	.042"/1.07mm	100cm	3.4cm	.038" (0.97mm)	RBR	bumper tip
---------	---	--------------	-------	-------	----------------	-----	------------

Amplatz Left Brachial/Radial


7703-B0	4	.042"/1.07mm	100cm	1.5cm	.038" (0.97mm)	AL2	bumper tip
---------	---	--------------	-------	-------	----------------	-----	------------



catalog number	French size	description		catheters and individual cat. no.
----------------	-------------	-------------	--	-----------------------------------

Cardiac Multipak

7751-10	4	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7704-10, 7700-20, 7701-20)	bumper tip
7751-20	4	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7704-20, 7700-20, 7701-20)	bumper tip
7560-13	5	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7529-11, 7521-21, 7523-21)	
7560-A3	5	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, PRC (7529-11, 7521-21, 7532-13)	
7560-23	5	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7529-23, 7521-21, 7523-21)	
7560-C3	5	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, PRC (7529-23, 7521-21, 7532-13)	
7560-33	5	CARDIAC MULTIPACK III	155° PIG (Braided), JL4.0, JR4.0 (7529-33, 7521-21, 7523-21)	
7560-43	5	CARDIAC MULTIPACK IV	PIG (Non-Braided), JL4.0, JR4.0 (7602-11, 7521-21, 7523-21)	
7560-53	5	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7602-23, 7521-21, 7523-21)	
7560-G3	5	CARDIAC MULTIPACK V	PIG (Braided), JL4.0 (high flow), JR4.0 (7529-11, 7601-21, 7523-21)	
7660-13	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7509-11, 7501-21, 7503-21)	
7660-A3	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, PRC (7509-11, 7501-21, 7512-13)	
7660-23	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7509-23, 7501-21, 7503-21)	
7660-B3	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, PRC (7509-23, 7501-21, 7512-13)	
7660-33	6	CARDIAC MULTIPACK III	155° PIG (Braided), JL4.0, JR4.0 (7509-33, 7501-21, 7503-21)	
7660-43	6	CARDIAC MULTIPACK IV	PIG (Non-Braided), JL4.0, JR4.0 (7569-11, 7501-21, 7503-21)	
7660-53	6	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7569-23, 7501-21, 7503-21)	
7660-63	6	CARDIAC MULTIPACK VI	155° PIG (Non-Braided), JL4.0, JR4.0 (7569-33, 7501-21, 7503-21)	



catalog number	French size	catheters (100cm)	sheath (11 cm) with .038 guide wire	guide wire (150cm) 3mm J .035" (0.89mm)Fixed Core
----------------	-------------	-------------------	-------------------------------------	---

Performa Multipak PLUS

4MPPRO-W-S	4	PRC, JL4.0, Straight Pigtail	4F PRO*	Yes
4MPPRO-W-145	4	JR4.0, JL4.0, 145° Pigtail	4F PRO*	Yes
5MPSW-S	5	JR 4.0, JL 4.0, Straight Pigtail	5F Prelude	Yes
5MPSW-145	5	JR 4.0, JL 4.0, 145° Pigtail	5F Prelude	Yes
5MPS-S	5	JR 4.0, JL 4.0, Straight Pigtail	5F Prelude	No
5MPS-145	5	JR 4.0, JL 4.0, 145° Pigtail	5F Prelude	No
5MPSW-SNB	5	JR 4.0, JL 4.0, Straight Pigtail (non braided)	5F Prelude	Yes
5MPPROW-S	5	JR 4.0, JL 4.0, Straight Pigtail	5F PRO*	Yes
5MPPRO-S	5	JR 4.0, JL 4.0, Straight Pigtail	5F PRO*	No
5MPPRO-145	5	JR 4.0, JL 4.0, 145° Pigtail	5F PRO*	No
5MPPRO-W-145	5	JR 4.0, JL 4.0, 145° Pigtail	5F PRO* **	Yes
5MPPROW-145	5	JR 4.0, JL 4.0, 145° Pigtail	5F PRO*	Yes
6MPSW-S	6	JR 4.0, JL 4.0, Straight Pigtail	6F Prelude	Yes
6MPSW-145	6	JR 4.0, JL 4.0, 145° PPigtail	6F Prelude	Yes
6MPS-S	6	JR 4.0, JL 4.0, Straight Pigtail	6F Prelude	No
6MPS-145	6	JR 4.0, JL 4.0, 145° Pigtail	6F Prelude	No
6MPPROW-S	6	JR 4.0, JL 4.0, Straight Pigtail	6F Pro*	Yes
6MPPROW-145	6	JR 4.0, JL 4.0, 145° Pigtail	6F Pro*	Yes









* Polypropylene (PRO) sheath tubing offers increased resilience and additional support during procedure.

** Contains 5F Prelude Pro™ with .035" guide wire.



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--

Pediatric Judkins Right and Left Catheters

7700-A0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 1.5	bumper tip	
7700-B0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 2.0	bumper tip	 JR 1.5
7700-C0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 2.5	bumper tip	
7700-D0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 3.0	bumper tip	 JR 2.0
7701-A0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 1.5	bumper tip	
7701-B0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 2.0	bumper tip	 JR 2.5
7701-C0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 2.5	bumper tip	
7701-D0	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 3.0	bumper tip	 JR 3.0
7521-L3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 1.5	bumper tip	 JL 1.5
7521-M3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 2.0	bumper tip	
7521-N3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 2.5	bumper tip	 JL 2.0
7521-P3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 3.0	bumper tip	 JL 2.5
7523-L3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 1.5	bumper tip	
7523-M3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 2.0	bumper tip	 JL 3.0
7523-N3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 2.5	bumper tip	
7523-P3	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 3.0	bumper tip	



For ordering information call Merit Medical customer service 1-800-35-MERIT.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT www.merit.com
Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40
Free phone for specific country:
Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71
Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

ID 051711 400463001/F