


Precision is TERUMO Territory.™

Runthrough® NS

Coronary Guidewires



Proven Terumo coronary guidewire technology

**Designed to be your first
choice coronary guidewire.**

 **TERUMO®**

TERUMO INTERVENTIONAL SYSTEMS

Precision you can feel.

Exceptional tip shape retention

- Aids vessel selection and lesion crossing
- Multi-vessel procedures with one wire
- Efficient, cost-effective procedures

DuoCore™ Technology




- 1:1 torque transfer
- Excellent steerability and pushability
- Fine tactile feel
- Excellent guide and device support

Optimal balance of hydrophobic and hydrophilic coatings

- Smooth trackability in tortuous vessels
- Superior device delivery

Tip durability¹


Tested to simulate guidewire experiencing resistance in a coronary vessel. Runthrough NS durable nitinol tip maintains its shape better than leading coronary guidewires.



1. Data on file. Terumo Corporation. Tokyo, Japan.



Runthrough[®] NS

Coronary Guidewires



Unique DuoCore Technology¹

Seamlessly fuses two proven guidewire technologies to ensure excellent steerability, pushability and trackability.



Balance MiddleWeight (BMW) is a registered trademark of Abbott Vascular. Prowater is a registered trademark of Asahi Intec Co.

TERUMO INTERVENTIONAL SYSTEMS

“The tip is worth its weight in gold...”

- Dr. Richard Heuser, Phoenix, AZ
Interventional Cardiologist

“The balloon feels as if it is floating down the wire with absolutely no resistance.”

- Dr. Lee Jobe, Raleigh, NC
Interventional Cardiologist

Terumo Runthrough® NS

Coronary Guidewires

Packed 5 per box. Re-shaper included.

Product Code	Diameter (in)	Length (cm)	Tip Radiopacity (cm)	Tip Shape
25-1011	0.014	180	3	Straight/Shapeable
25-1013	0.014	300	3	Straight/Shapeable

Packed 1 per box. Does not include a re-shaper.

25-9092	0.014	150	–	Extension wire
---------	-------	-----	---	----------------

TERUMO—A global leader in medical devices

Terumo Sheaths
Terumo Guidewires
Terumo Angiographic Catheters
Terumo Coronary Catheters
Terumo Microcatheters
Interventional Products

For more information, call 800.862.4143 to speak to an Inside Sales Customer Care Representative, or visit www.terumo/is.com

Terumo Interventional Systems · 2101 Cottontail Lane · Somerset, NJ 08873 · Fax: 800.411.5870

©2009 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. Lit. # TIS100r1 07/09